

For Research Use Only

MATN3 Recombinant antibody, PBS Only

Catalog Number: 86125-3-PBS



Basic Information

Catalog Number:

86125-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3136

GenBank Accession Number:

NM_002381.4

GeneID (NCBI):

4148

UNIPROT ID:

O15232

Full Name:

matrilin 3

Calculated MW:

53 kDa

Observed MW:

53 kDa

Purification Method:

Protein A purification

CloneNo.:

250692E1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Matrilin-3 is a protein that in humans is encoded by the MATN3 gene. It is linked to the development of many types of cartilage, and part of the Matrilin family, which includes Matrilin-1, Matrilin-2, Matrilin-3, and Matrilin-4, a family of filamentous-forming adapter oligomeric extracellular proteins that are linked to the formation of cartilage and bone, as well as maintaining homeostasis after development. It is considered an extracellular matrix protein that functions as an adapter protein where the Matrilin-3 subunit can form both homo-tetramers and hetero-oligomers with subunits from Matrilin-1 which is the cartilage matrix protein. This restricted tissue has been strongly expressed in growing skeletal tissue as well as cartilage and bone. (PMID:26499313, 27104523, 11950247)

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86125-3-RR (MATN3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86125-3-PBS in a different storage buffer formulation.