

For Research Use Only

NRP2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86126-1-PBS



Basic Information

Catalog Number:

86126-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3125

GenBank Accession Number:

NM_201266.1

GeneID (NCBI):

8828

UNIPROT ID:

O60462-1

Full Name:

neuropilin 2

Calculated MW:

105 kDa

Purification Method:

Protein A purification

CloneNo.:

250688E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Neuropilin 2 concentration was determined to be 15.7 ng/mL in human serum.



Cytometric bead array standard curve of MP02281-1, Neuropilin 2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86126-2-PBS. Detection antibody: 86126-1-PBS. Standard: Eg3125. Range: 0.781-100 ng/mL.