

For Research Use Only

# PUM1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86129-2-PBS



## Basic Information

Catalog Number:

86129-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24651

GenBank Accession Number:

BC013398

GeneID (NCBI):

9698

UNIPROT ID:

Q14671

Full Name:

pumilio homolog 1 (Drosophila)

Calculated MW:

1186 aa, 126 kDa

Purification Method:

Protein A purification

CloneNo.:

250792F11

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean PUM1 concentration was determined to be 2.1 ng/mL in MCF-7 cell extract based on a 1.3 mg/mL extract load.



Sandwich ELISA standard curve of MP02520-1, Human PUM1 Recombinant Matched Antibody Pair - PBS only. 86129-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag24651. 86129-3-PBS was HRP conjugated as the detection antibody. Range: 156-5000 pg/mL