

For Research Use Only

PUM1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86129-3-PBS



Basic Information

Catalog Number: 86129-3-PBS	GenBank Accession Number: BC013398	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 9698	CloneNo.: 250792E8
Isotype: IgG	UNIPROT ID: Q14671	
Immunogen Catalog Number: AG24651	Full Name: pumilio homolog 1 (Drosophila)	
	Calculated MW: 1186 aa, 126 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PUM1 concentration was determined to be 2.1 ng/mL in MCF-7 cell extract based on a 1.3 mg/mL extract load.



Sandwich ELISA standard curve of MP02520-1, Human PUM1 Recombinant Matched Antibody Pair - PBS only. 86129-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag24651. 86129-3-PBS was HRP conjugated as the detection antibody. Range: 156-5000 pg/mL