

For Research Use Only

GMPS Recombinant antibody

Catalog Number: 86130-1-RR



Basic Information

Catalog Number:

86130-1-RR

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9350

GenBank Accession Number:

BC012178

GeneID (NCBI):

8833

UNIPROT ID:

P49915

Full Name:

guanine monphosphate synthetase

Calculated MW:

693 aa, 77 kDa

Observed MW:

67 kDa

Purification Method:

Protein A purification

CloneNo.:

250704D4

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, Jurkat cells, HL-60 cells, Ramos cells, mouse testis tissue, rat testis tissue, mouse brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86130-1-RR (GMPS antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.