

For Research Use Only

GMPS Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86130-2-PBS



Basic Information

Catalog Number:

86130-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9350

GenBank Accession Number:

BC012178

GeneID (NCBI):

8833

UNIPROT ID:

P49915

Full Name:

guanine monphosphate synthetase

Calculated MW:

693 aa, 77 kDa

Purification Method:

Protein A purification

CloneNo.:

250704B4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean GMPS concentration was determined to be 12.2 ng/mL in HeLa cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP02408-1, GMPS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86130-3-PBS. Detection antibody: 86130-2-PBS. Standard: Ag9350. Range: 0.469-60 ng/mL.