

For Research Use Only

BTNL2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86137-3-PBS



Basic Information

Catalog Number: 86137-3-PBS	GenBank Accession Number: BC127642	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 56244	CloneNo.: 250562C10
Source: Rabbit	UNIPROT ID: Q9UIRO	
Isotype: IgG	Full Name: butyrophilin-like 2 (MHC class II associated)	
Immunogen Catalog Number: AG18803	Calculated MW: 455 aa, 50 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02518-1, Human BTNL2 Recombinant Matched Antibody Pair - PBS only. 86137-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag18803. 86137-3-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



The mean BTNL2 concentration was determined to be 39.10 ng/mL in HeLa cell extract based on a 1.1 mg/mL extract load.