

For Research Use Only

# DAB2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86141-3-PBS



## Basic Information

Catalog Number:

86141-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0162

GenBank Accession Number:

BC003064

GeneID (NCBI):

1601

UNIPROT ID:

P98082

Full Name:

disabled homolog 2, mitogen-responsive phosphoprotein (Drosophila)

Calculated MW:

82 kDa

Purification Method:

Protein A purification

CloneNo.:

250409C12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean DAB2 concentration was determined to be 34.8 ng/mL in Hela cell extract based on a 1.1 mg/mL extract load.



Cytometric bead array standard curve of MP02525-1, DAB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86141-1-PBS. Detection antibody: 86141-3-PBS. Standard: Ag0162. Range: 0.625-80 ng/mL.