

For Research Use Only

LILRB4/CD85k/ILT3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86142-2-PBS



Basic Information

Catalog Number:

86142-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1743

GenBank Accession Number:

AAH26309.1

GeneID (NCBI):

11006

UNIPROT ID:

Q8NHJ6

Full Name:

leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 4

Calculated MW:

49kDa

Purification Method:

Protein A purification

CloneNo.:

250356F6

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02395-1, Human LILRB4/CD85k/ILT3 Recombinant Matched Antibody Pair - PBS only. 86142-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1743. 86142-2-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.



The mean LILRB4/CD85k/ILT3 concentration was determined to be 6,086.5 pg/mL in A549 cell extract based on a 1.2 mg/mL extract load.