

For Research Use Only

AKR1A1 Recombinant antibody

Catalog Number: 86147-1-RR



Basic Information

Catalog Number:

86147-1-RR

Concentration:

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7037

GenBank Accession Number:

BC000670

GeneID (NCBI):

10327

UNIPROT ID:

P14550

Full Name:

aldo-keto reductase family 1,
member A1 (aldehyde reductase)

Calculated MW:

37 kDa

Observed MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

250713E8

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, mouse lung tissue, HepG2 cells, L02
cells, HCT 116 cells, LNCaP cells, rat lung tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86147-1-RR (AKR1A1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 86147-1-RR (AKR1A1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.