

For Research Use Only

TP63 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86154-4-PBS



Basic Information

Catalog Number:

86154-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2791

GenBank Accession Number:

BC039815

GeneID (NCBI):

8626

UNIPROT ID:

Q9H3D4

Full Name:

tumor protein p63

Calculated MW:

680 aa, 77 kDa

Purification Method:

Protein A purification

CloneNo.:

250042A4

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02275-1, Human p63 Recombinant Matched Antibody Pair - PBS only. 86154-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag2791. 86154-3-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL



The mean p63 concentration was determined to be 236.27 ng/mL in A431 cell extract based on a 3.7 mg/mL extract load and 25.49 ng/mL in HaCat cell extract based on a 1.2 mg/mL extract load.