

For Research Use Only

MYH10 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86157-2-PBS



Basic Information

Catalog Number:

86157-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16113

GenBank Accession Number:

BC150634

GeneID (NCBI):

4628

UNIPROT ID:

P35580

Full Name:

myosin, heavy chain 10, non-muscle

Calculated MW:

1985 aa, 230 kDa

Purification Method:

Protein A purification

CloneNo.:

250791E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MYH10 concentration was determined to be 22.2 ng/mL in HepG2 cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02284-1, MYH10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86157-2-PBS. Detection antibody: 86157-1-PBS. Standard: Ag16113. Range: 0.625-80 ng/mL.