

For Research Use Only

Mouse Mcam Recombinant antibody, PBS Only (Detector)

Catalog Number: 86159-1-PBS



Basic Information

Catalog Number:

86159-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2667

GenBank Accession Number:

NM_023061.2

GeneID (NCBI):

84004

UNIPROT ID:

Q8R2Y2-1

Full Name:

melanoma cell adhesion molecule

Calculated MW:

72 kDa

Purification Method:

Protein A purification

CloneNo.:

250806F6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD146/MCAM concentration was determined to be 2035.9 ng/mL in mouse lung tissue extract based on a 5.9 mg/mL extract load.



Cytometric bead array standard curve of MP02286-1, MOUSE CD146/MCAM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86159-2-PBS. Detection antibody: 86159-1-PBS. Standard: Eg2667. Range: 0.781-100 ng/mL