

For Research Use Only

REEP4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86164-5-PBS



Basic Information

Catalog Number:

86164-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24569

GenBank Accession Number:

BC013048

GeneID (NCBI):

80346

UNIPROT ID:

Q9H6H4

Full Name:

receptor accessory protein 4

Calculated MW:

29 kDa

Purification Method:

Protein A purification

CloneNo.:

250770D9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean REEP4 concentration was determined to be 5.4 ng/mL in HeLa cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP02417-1, REEP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86164-5-PBS. Detection antibody: 86164-4-PBS. Standard: Ag24569. Range: 0.781-100 ng/mL.