

For Research Use Only

# PLEKHH2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86166-2-PBS



## Basic Information

Catalog Number:

86166-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5409

GenBank Accession Number:

BC063310

GeneID (NCBI):

130271

UNIPROT ID:

Q8IVE3

Full Name:

pleckstrin homology domain  
containing, family H (with MyTH4  
domain) member 2

Calculated MW:

168 kDa

Purification Method:

Protein A purification

CloneNo.:

250291C2

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP02268-1, Human PLEKHH2 Recombinant Matched Antibody Pair - PBS only. 86166-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag5409. 86166-1-PBS was HRP conjugated as the detection antibody. Range: 0.625-20 ng/mL



The mean PLEKHH2 concentration was determined to be 5.99 ng/mL in Jurkat cell extract based on a 1.90 mg/mL extract load.