

For Research Use Only

CSNK2A2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86169-3-PBS



Basic Information

Catalog Number: 86169-3-PBS	GenBank Accession Number: BC008812	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 1459	CloneNo.: 250812F8
Isotype: IgG	UNIPROT ID: P19784	
Immunogen Catalog Number: AG0877	Full Name: casein kinase 2, alpha prime polypeptide	
	Calculated MW: 41 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CSNK2A2 concentration was determined to be 23.9 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02423-1, CSNK2A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86169-4-PBS. Detection antibody: 86169-3-PBS. Standard: Ag0877. Range: 0.469-60 ng/mL.