

For Research Use Only

# CDK5RAP3 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86177-3-PBS



## Basic Information

Catalog Number:

86177-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1474

GenBank Accession Number:

BC009957

GeneID (NCBI):

80279

UNIPROT ID:

Q96JB5

Full Name:

CDK5 regulatory subunit associated protein 3

Calculated MW:

57 kDa

Purification Method:

Protein A purification

CloneNo.:

250801A6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean CDK5RAP3 concentration was determined to be 122.4 ng/mL in A431 cell extract based on a 3.7 mg/mL extract load.



Cytometric bead array standard curve of MP02384-1, CDK5RAP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86177-3-PBS. Detection antibody: 86177-1-PBS. Standard: Ag1474. Range: 0.469-60 ng/mL.