

For Research Use Only

BCL2L13 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86183-3-PBS



Basic Information

Catalog Number: 86183-3-PBS	GenBank Accession Number: BC003032	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 23786	CloneNo.: 250705F8
Isotype: IgG	UNIPROT ID: Q9BXK5	
Immunogen Catalog Number: AG9841	Full Name: BCL2-like 13 (apoptosis facilitator)	
	Calculated MW: 485 aa, 53 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean BCL2L13 concentration was determined to be 2.9 ng/mL in HepG2 cell extract based on a 1.1 mg/mL extract load.



Cytometric bead array standard curve of MP02404-1, BCL2L13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86183-3-PBS. Detection antibody: 86183-2-PBS. Standard: Ag9841. Range: 0.781-100 ng/mL.



Bi-layer interferometry (BLI) kinetic assay of 86183-3-PBS against Human BCL2L13 was performed. The affinity constant is 1.27 nM.