

For Research Use Only

GRAMD1B Recombinant antibody, PBS Only (Detector)

Catalog Number: 86185-2-PBS



Basic Information

Catalog Number:

86185-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20670

GenBank Accession Number:

BC107480

GeneID (NCBI):

57476

UNIPROT ID:

Q3KR37

Full Name:

GRAM domain containing 1B

Calculated MW:

738 aa, 85 kDa

Purification Method:

Protein A purification

CloneNo.:

250789B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean GRAMD1B concentration was determined to be 2.0 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02407-1, GRAMD1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86185-1-PBS. Detection antibody: 86185-2-PBS. Standard: Ag20670. Range: 0.469-60 ng/mL.