

For Research Use Only

CXADR Recombinant antibody, PBS Only

Catalog Number: 86186-1-PBS



Basic Information

Catalog Number:

86186-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34346

GenBank Accession Number:

BC003684

GeneID (NCBI):

1525

UNIPROT ID:

P78310

Full Name:

coxsackie virus and adenovirus
receptor

Calculated MW:

365 aa, 40 kDa

Observed MW:

45-55 kDa

Purification Method:

Protein A purification

CloneNo.:

250811C8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86186-1-RR (CXADR antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86186-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 86186-1-RR (CXADR antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86186-1-PBS in a different storage buffer formulation.