

For Research Use Only

CXADR Recombinant antibody

Catalog Number: 86186-1-RR



Basic Information

Catalog Number:

86186-1-RR

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34346

GenBank Accession Number:

BC003684

GeneID (NCBI):

1525

UNIPROT ID:

P78310

Full Name:

coxsackie virus and adenovirus
receptor

Calculated MW:

365 aa, 40 kDa

Observed MW:

45-55 kDa

Purification Method:

Protein A purification

CloneNo.:

250811C8

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HaCaT cells, mouse heart tissue, rat heart tissue,
mouse small intestine tissue, BxPC-3 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86186-1-RR (CXADR antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 86186-1-RR (CXADR antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.