

For Research Use Only

CXADR Recombinant antibody, PBS Only (Detector)

Catalog Number: 86186-2-PBS



Basic Information

Catalog Number:

86186-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34346

GenBank Accession Number:

BC003684

GeneID (NCBI):

1525

UNIPROT ID:

P78310

Full Name:

coxsackie virus and adenovirus
receptor

Calculated MW:

365 aa, 40 kDa

Purification Method:

Protein A purification

CloneNo.:

250811C7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CXADR concentration was determined to be 60.4 ng/mL in HeLa cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02296-1, CXADR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86186-3-PBS. Detection antibody: 86186-2-PBS. Standard: Ag34346. Range: 1.563-200 ng/mL.