

For Research Use Only

# CXADR Recombinant antibody, PBS Only (Capture)

Catalog Number: 86186-3-PBS



## Basic Information

Catalog Number:

86186-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34346

GenBank Accession Number:

BC003684

GeneID (NCBI):

1525

UNIPROT ID:

P78310

Full Name:

coxsackie virus and adenovirus  
receptor

Calculated MW:

365 aa, 40 kDa

Purification Method:

Protein A purification

CloneNo.:

250811F9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean CXADR concentration was determined to be 60.4 ng/mL in HeLa cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02296-1, CXADR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86186-3-PBS. Detection antibody: 86186-2-PBS. Standard: Ag34346. Range: 1.563-200 ng/mL.