

For Research Use Only

CXADR Recombinant antibody, PBS Only (Detector)

Catalog Number: 86186-4-PBS



Basic Information

Catalog Number:

86186-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34346

GenBank Accession Number:

BC003684

GeneID (NCBI):

1525

UNIPROT ID:

P78310

Full Name:

coxsackie virus and adenovirus
receptor

Calculated MW:

365 aa, 40 kDa

Purification Method:

Protein A purification

CloneNo.:

250811D9

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02296-2, Human CXADR Recombinant Matched Antibody Pair - PBS only. 86186-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag34346. 86186-4-PBS was HRP conjugated as the detection antibody. Range: 0.234-15 ng/mL.



The mean CXADR concentration was determined to be 26.27 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load and 183.10 ng/mL in BxPC-3 cell extract based on a 1.20 mg/mL extract load.