

For Research Use Only

FKBP5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86189-3-PBS



Basic Information

Catalog Number:

86189-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5337

GenBank Accession Number:

BC042605

GeneID (NCBI):

2289

UNIPROT ID:

Q13451

Full Name:

FK506 binding protein 5

Calculated MW:

457 aa, 51 kDa

Purification Method:

Protein A purification

CloneNo.:

250798C11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FKBP5 concentration was determined to be 473.3 ng/mL in K-562 cell extract based on a 3.00 mg/mL extract load.



Cytometric bead array standard curve of MP02289-1, FKBP5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86189-3-PBS. Detection antibody: 86189-2-PBS. Standard: Ag5337. Range: 1.563-200 ng/mL.