

For Research Use Only

TNK1 Recombinant antibody, PBS Only

Catalog Number: 86190-2-PBS



Basic Information

Catalog Number:

86190-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5402

GenBank Accession Number:

BC035782

GeneID (NCBI):

8711

UNIPROT ID:

Q13470

Full Name:

tyrosine kinase, non-receptor, 1

Calculated MW:

72 kDa

Observed MW:

72 kDa

Purification Method:

Protein A purification

CloneNo.:

250786A6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86190-2-RR (TNK1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86190-2-PBS in a different storage buffer formulation.



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 86190-2-RR (TNK1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86190-2-PBS in a different storage buffer formulation.