

For Research Use Only

TNK1 Recombinant antibody

Catalog Number: 86190-2-RR



Basic Information

Catalog Number:

86190-2-RR

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5402

GenBank Accession Number:

BC035782

GeneID (NCBI):

8711

UNIPROT ID:

Q13470

Full Name:

tyrosine kinase, non-receptor, 1

Calculated MW:

72 kDa

Observed MW:

72 kDa

Purification Method:

Protein A purification

CloneNo.:

250786A6

Recommended Dilutions:

WB: 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : K-562 cells, fetal human brain tissue, PC-3 cells, HaCaT cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86190-2-RR (TNK1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 86190-2-RR (TNK1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.