

For Research Use Only

RUNX3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86191-3-PBS



Basic Information

Catalog Number: 86191-3-PBS	GenBank Accession Number: BC013362	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 864	CloneNo.: 250589D4
Source: Rabbit	UNIPROT ID: Q13761	
Isotype: IgG	Full Name: runt-related transcription factor 3	
Immunogen Catalog Number: AG24339	Calculated MW: 44 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02431-1, Human RUNX3 Recombinant Matched Antibody Pair - PBS only. 86191-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag24339. 86191-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL

The mean RUNX3 concentration was determined to be 76.82 ng/mL in Jurkat cell extract based on a 1.20 mg/mL extract load.