

For Research Use Only

ANGPTL3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86193-2-PBS



Basic Information

Catalog Number:

86193-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2573

GenBank Accession Number:

BC058287

GeneID (NCBI):

27329

UNIPROT ID:

Q9Y5C1

Full Name:

angiopoietin-like 3

Calculated MW:

460 aa, 54 kDa

Purification Method:

Protein A purification

CloneNo.:

250808A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ANGPTL3 concentration was determined to be 14.8 ng/mL in human serum.



Cytometric bead array standard curve of MP02294-1, ANGPTL3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86193-2-PBS. Detection antibody: 86193-1-PBS. Standard: Ag2573. Range: 1.563-200 ng/mL.