

For Research Use Only

OLR1/LOX1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86194-2-PBS



Basic Information

Catalog Number:	86194-2-PBS	GenBank Accession Number:	NM_002543.3	Purification Method:	Protein A purification
Concentration:	1 mg/ml	GeneID (NCBI):	4973	CloneNo.:	250525B12
Source:	Rabbit	UNIPROT ID:	P78380-1		
Isotype:	IgG	Full Name:	oxidized low density lipoprotein (lectin-like) receptor 1		
Immunogen Catalog Number:	EG2911	Calculated MW:	31 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean OLR1/LOX1 concentration was determined to be 5.0 ng/mL in HeLa cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02293-1, OLR1/LOX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86194-2-PBS. Detection antibody: 86194-1-PBS. Standard: Eg2911. Range: 0.156-20 ng/mL.