

For Research Use Only

Frataxin Recombinant antibody, PBS Only (Detector)

Catalog Number: 86195-1-PBS



Basic Information

Catalog Number:

86195-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4622

GenBank Accession Number:

NM_000144.5

GeneID (NCBI):

2395

UNIPROT ID:

Q16595-1

Full Name:

frataxin

Calculated MW:

23kDa

Purification Method:

Protein A purification

CloneNo.:

250761A7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Frataxin concentration was determined to be 779.5 ng/mL in HepG2 cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02292-1, Frataxin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86195-2-PBS. Detection antibody: 86195-1-PBS. Standard: Eg4622. Range: 0.313-40 ng/mL.