

For Research Use Only

CBFA2T3 Recombinant antibody, PBS Only

Catalog Number: 86198-4-PBS



Basic Information

Catalog Number:

86198-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9265

GenBank Accession Number:

BC062624

GeneID (NCBI):

863

UNIPROT ID:

O75081

Full Name:

core-binding factor, runt domain,
alpha subunit 2; translocated to, 3

Calculated MW:

71 kDa

Observed MW:

71 kDa

Purification Method:

Protein A purification

CloneNo.:

250859C5

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86198-4-RR (CBFA2T3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86198-4-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed Jurkat cells using CBFA2T3 antibody (86198-4-RR, Clone: 250859C5) at dilution of 1:500 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 86198-4-PBS in a different storage buffer formulation.