

For Research Use Only

# DUSP9 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86199-1-PBS



## Basic Information

Catalog Number: 86199-1-PBS	GenBank Accession Number: BC060837	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 1852	CloneNo.: 250824E8
Source: Rabbit	UNIPROT ID: Q99956	
Isotype: IgG	Full Name: dual specificity phosphatase 9	
Immunogen Catalog Number: AG25167	Calculated MW: 42 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



The mean DUSP9 concentration was determined to be 4.0 ng/mL in HepG2 cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02300-1, DUSP9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86199-2-PBS. Detection antibody: 86199-1-PBS. Standard: Ag25167. Range: 0.156-20 ng/mL.