

For Research Use Only

DUSP9 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86199-2-PBS



Basic Information

Catalog Number: 86199-2-PBS	GenBank Accession Number: BC060837	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 1852	CloneNo.: 250824E12
Source: Rabbit	UNIPROT ID: Q99956	
Isotype: IgG	Full Name: dual specificity phosphatase 9	
Immunogen Catalog Number: AG25167	Calculated MW: 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean DUSP9 concentration was determined to be 4.0 ng/mL in HepG2 cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02300-1, DUSP9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86199-2-PBS. Detection antibody: 86199-1-PBS. Standard: Ag25167. Range: 0.156-20 ng/mL.