

For Research Use Only

DUSP9 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86199-4-PBS



Basic Information

Catalog Number:

86199-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25167

GenBank Accession Number:

BC060837

GeneID (NCBI):

1852

UNIPROT ID:

Q99956

Full Name:

dual specificity phosphatase 9

Calculated MW:

42 kDa

Purification Method:

Protein A purification

CloneNo.:

250824B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean DUSP9 concentration was determined to be 1.9 ng/mL in HepG2 cell extract based on a 1.0 mg/mL extract load, Raji cell extract were served as a negative control.



Cytometric bead array standard curve of MP02300-2, DUSP9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86199-5-PBS. Detection antibody: 86199-4-PBS. Standard: Ag25167. Range: 0.156-20 ng/mL.