

For Research Use Only

PSMD12 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86200-2-PBS



Basic Information

Catalog Number: 86200-2-PBS	GenBank Accession Number: BC019062	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 5718	CloneNo.: 250848G4
Source: Rabbit	UNIPROT ID: O00232	
Isotype: IgG	Full Name: proteasome (prosome, macropain) 26S subunit, non-ATPase, 12	
Immunogen Catalog Number: AG1957	Calculated MW: 456 aa, 53 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PSMD12 concentration was determined to be 45.7 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02301-1, PSMD12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86200-2-PBS. Detection antibody: 86200-1-PBS. Standard: Ag1957. Range: 1.563-200 ng/mL.