

For Research Use Only

Mouse Dkk-1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86207-2-PBS



Basic Information

Catalog Number:

86207-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4649

GenBank Accession Number:

NM_010051.3

GeneID (NCBI):

13380

UNIPROT ID:

O54908

Full Name:

dickkopf homolog 1 (Xenopus laevis)

Calculated MW:

29kDa

Purification Method:

Protein A purification

CloneNo.:

250658A3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02491-1, Mouse DKK1 Recombinant Matched Antibody Pair - PBS only. 86207-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg4649. 86207-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



Serum of eight mice was measured. The mouse DKK1 concentration of detected samples was determined to be 4.72 ng/mL with a range of 2.40 - 7.06 ng/mL