

For Research Use Only

STK25 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86216-2-PBS



Basic Information

Catalog Number:

86216-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22964

GenBank Accession Number:

BC015793

GeneID (NCBI):

10494

UNIPROT ID:

O00506

Full Name:

serine/threonine kinase 25 (STE20 homolog, yeast)

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

250202G5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02398-1, Human STK25 Recombinant Matched Antibody Pair - PBS only. 86216-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag22964. 86216-1-PBS was HRP conjugated as the detection antibody. Range: 62.5-4000 pg/mL



The mean SNAI2/SLUG concentration was determined to be 4,108.0 pg/mL in LnCap cell extract based on a 2.0 mg/mL extract load and 3,191.3 pg/mL in U2OS cell extract based on a 2.2 mg/mL extract load.