

For Research Use Only

HOXC10 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86219-2-PBS



Basic Information

Catalog Number:

86219-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2655

GenBank Accession Number:

BC001293

GeneID (NCBI):

3226

UNIPROT ID:

Q9NYD6

Full Name:

homeobox C10

Calculated MW:

342 aa, 38 kDa

Purification Method:

Protein A purification

CloneNo.:

250873D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean HOXC10 concentration was determined to be 5.3 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02297-1, HOXC10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86219-2-PBS. Detection antibody: 86219-1-PBS. Standard: Ag2655. Range: 0.781-100 ng/mL.