

For Research Use Only

# TORC1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86226-2-PBS



## Basic Information

Catalog Number:

86226-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0693

GenBank Accession Number:

BC028050

GeneID (NCBI):

23373

UNIPROT ID:

Q6UUV9

Full Name:

CREB regulated transcription  
coactivator 1

Calculated MW:

75 kDa

Purification Method:

Protein A purification

CloneNo.:

250870F5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean TORC1/CRTC1 concentration was determined to be 209.6 ng/mL in Jurkat cell extract based on a 2.2 mg/mL extract load.



Cytometric bead array standard curve of MP02309-1, TORC1/CRTC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86226-1-PBS. Detection antibody: 86226-2-PBS. Standard: Ag0693. Range: 0.781-100 ng/mL.