

For Research Use Only

CUGBP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86227-2-PBS



Basic Information

Catalog Number:

86227-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3677

GenBank Accession Number:

BC031079

GeneID (NCBI):

10658

UNIPROT ID:

Q92879

Full Name:

CUG triplet repeat, RNA binding protein 1

Calculated MW:

483 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

250200B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CUGBP1 concentration was determined to be 7.3 ng/mL in HL-60 cell extract based on a 3.7 mg/mL extract load.



Cytometric bead array standard curve of MP02347-1, CUGBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86227-2-PBS. Detection antibody: 86227-1-PBS. Standard: Ag3677. Range: 1.562-100 ng/mL.