

For Research Use Only

POLDIP2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86232-1-PBS



Basic Information

Catalog Number:

86232-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7154

GenBank Accession Number:

BC000655

GeneID (NCBI):

26073

UNIPROT ID:

Q9Y257

Full Name:

polymerase (DNA-directed), delta interacting protein 2

Calculated MW:

42 kDa

Purification Method:

Protein A purification

CloneNo.:

250999F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean POLDIP2 concentration was determined to be 84.6 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02345-1, POLDIP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86232-2-PBS. Detection antibody: 86232-1-PBS. Standard: Ag7154. Range: 0.781-100 ng/mL.