

For Research Use Only

PSMC1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86237-3-PBS



Basic Information

Catalog Number: 86237-3-PBS	GenBank Accession Number: BC016368	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 5700	CloneNo.: 251014F11
Source: Rabbit	UNIPROT ID: P62191	
Isotype: IgG	Full Name: proteasome (prosome, macropain) 26S subunit, ATPase, 1	
Immunogen Catalog Number: AG1677	Calculated MW: 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PSMC1 concentration was determined to be 1.8 ng/mL in HeLa cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02311-1, PSMC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86237-3-PBS. Detection antibody: 86237-2-PBS. Standard: Ag1677. Range: 0.156-20 ng/mL.