

For Research Use Only

Pan-Keratin (Type I) Recombinant monoclonal antibody, PBS Only

Catalog Number: 86238-2-PBS



Basic Information

Catalog Number:

86238-2-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Pan-Keratin (Type I)

UNIPROT ID:

Q15323

Full Name:

Observed MW:

45-65 kDa

Purification Method:

Protein A purification

CloneNo.:

251039C1

Applications

Tested Applications:

WB, IHC, IF-P, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Keratin is a family of fibrous structural proteins that are key components in hair, nails, feathers, hooves, and the outer layer of skin (epidermis). There are two main types of keratin: Type I (acidic) and Type II (basic/neutral), which form heterodimers to create intermediate filaments, providing mechanical strength to cells. Type I keratins (or Type I cytokeratins) are acidic, low molecular weight (40-64 kDa) proteins that form Type I intermediate filaments (IFs) in the cytoplasmic cytoskeleton of all mammalian epithelial cells. They are encoded on chromosome 17q and include keratins K9-K10, K12-K28, and K31-K40.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 86238-2-RR (Pan-Keratin (Type I) antibody) at dilution of 1:30000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86238-2-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 86238-2-RR (Pan-Keratin (Type I) antibody) at dilution of 1:1000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 86238-2-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 86238-2-RR (Pan-Keratin (Type I) antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 86238-2-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded human cervical cancer tissue using Pan-Keratin (Type I) antibody (86238-2-RR, Clone: 251039C1) at dilution of 1:500 and CoraLite® 594-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-4), Pan-Keratin (Type II) antibody 86370-1-RR labeled with FlexAble 2.0 CoraLite® Plus 488 Kit (KFA501, green). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using

