

For Research Use Only

PDIA4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86239-3-PBS



Basic Information

Catalog Number: 86239-3-PBS	GenBank Accession Number: NM_004911.5	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 9601	CloneNo.: 250697H3
Source: Rabbit	UNIPROT ID: P13667	
Isotype: IgG	Full Name: protein disulfide isomerase family A, member 4	
Immunogen Catalog Number: EG3011	Calculated MW: 73 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ERp72 concentration was determined to be 3418.9 ng/mL in A431 cell extract based on a 3.7 mg/mL extract load.



Cytometric bead array standard curve of MP02383-1, ERp72 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86239-3-PBS. Detection antibody: 86239-2-PBS. Standard: Eg3011. Range: 0.625-40 ng/mL.