

For Research Use Only

Draxin Recombinant antibody, PBS Only (Capture)

Catalog Number: 86244-2-PBS



Basic Information

Catalog Number: 86244-2-PBS	GenBank Accession Number: NM_198545.3	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 374946	CloneNo.: 250677C11
Source: Rabbit	UNIPROT ID: Q8NBI3	
Isotype: IgG	Full Name: chromosome 1 open reading frame 187	
Immunogen Catalog Number: EG4203	Calculated MW: 39 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02569-1, Human Draxin Recombinant Matched Antibody Pair - PBS only. 86244-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg4203. 86244-3-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL.



The mean Draxin concentration was determined to be 30.83 ng/mL in SH-SY5Y cell extract based on a 1.40 mg/mL extract load.