

For Research Use Only

Rat CD36/SR-B3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86246-3-PBS



Basic Information

Catalog Number:

86246-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3222

GenBank Accession Number:

GeneID (NCBI):

29184

UNIPROT ID:

Q07969

Full Name:

CD36 molecule (thrombospondin receptor)

Calculated MW:

53 kDa

Purification Method:

Protein A purification

CloneNo.:

250237D5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD36 concentration was determined to be 2888.5 ng/mL in Rat spleen tissue extract based on a 8.8 mg/mL extract load.



Cytometric bead array standard curve of MP02307-1, RAT CD36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86246-3-PBS. Detection antibody: 86246-2-PBS. Standard: Eg3222. Range: 0.313-40 ng/mL.