

For Research Use Only

BAG5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86248-2-PBS



Basic Information

Catalog Number:	86248-2-PBS	GenBank Accession Number:	BC044216	Purification Method:	Protein A purification
Concentration:	1 mg/ml	GeneID (NCBI):	9529	CloneNo.:	250877G7
Source:	Rabbit	UNIPROT ID:	Q9UL15		
Isotype:	IgG	Full Name:	BCL2-associated athanogene 5		
Immunogen Catalog Number:	AG24320	Calculated MW:	51 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BAG5 concentration was determined to be 110.8 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02328-1, BAG5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86248-3-PBS. Detection antibody: 86248-2-PBS. Standard: Ag24320. Range: 0.781-100 ng/mL.