

For Research Use Only

FKBP1A Recombinant antibody, PBS Only (Capture)

Catalog Number: 86249-3-PBS



Basic Information

Catalog Number:

86249-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0406

GenBank Accession Number:

BC001925

GeneID (NCBI):

2280

UNIPROT ID:

P62942

Full Name:

FK506 binding protein 1A, 12kDa

Calculated MW:

12 kDa

Purification Method:

Protein A purification

CloneNo.:

250923D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FKBP1A concentration was determined to be 149.2 ng/mL in Jurkat cell extract based on a 2.2 mg/mL extract load.



Cytometric bead array standard curve of MP02313-1, FKBP1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86249-3-PBS. Detection antibody: 86249-2-PBS. Standard: Ag0406. Range: 0.781-100 ng/mL.