

For Research Use Only

# KRT19 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86250-7-PBS



## Basic Information

Catalog Number:

86250-7-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC002539

GeneID (NCBI):

3880

UNIPROT ID:

P08727

Full Name:

keratin 19

Calculated MW:

44 kDa

Purification Method:

Protein A purification

CloneNo.:

250894F4

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP02308-2, Human Cytokeratin 19 Recombinant Matched Antibody Pair - PBS only. 86250-8-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY02555. 86250-7-PBS was HRP conjugated as the detection antibody. Range: 1.56-100 ng/mL



The mean Cytokeratin 19 concentration was determined to be 228.05 ng/mL in MCF-7 cell extract based on a 1.30 mg/mL extract load.